

You can make a difference!

After flossing and brushing treat your teeth with GC Tooth Mousse and help everyone in the family to achieve a healthy smile.

GC Tooth Mousse contains 10% (wt) RECALDENT™ (CPP-ACP) to provide you with a concentrated source of bio-available calcium and phosphate for strengthening and protection of teeth. And it tastes delicious!

So easy to use!

1



Squeeze a small amount of GC Tooth Mousse onto your finger.

2



Apply onto your teeth and use your tongue to spread evenly around all teeth.

3

For maximum benefit, leave GC Tooth Mousse on the teeth for as long as possible. The minimum recommended application time is three minutes. At the end of the application you can expectorate, rinse or swallow the remainder.

GC Tooth Mousse comes in individual 40g tubes in a variety of delicious flavours.



GC

Mi

Recaldent

GC Asia Dental Pte. Ltd
Changi Logistics Centre
19 Loyang Way #06-27 Singapore 508724
Tel + 65 6546 7588 Fax + 65 6546 7577
www.gcasia.info © GC Corporation 2010

GC Tooth Mousse contains RECALDENT™ (CPP-ACP), a unique ingredient developed at The School of Dental Science, The University of Melbourne, Victoria, Australia. RECALDENT and RECALDENT Device are trademarks used under licence. GC Tooth Mousse should not be used by people with milk protein allergies. If any allergic reaction occurs, this may indicate sensitivity to the benzoate preservatives, or to some other component of the product. In this event, discontinue use of the product and contact your physician.

HS Code: TM-BROCHILD

GC

GC Tooth Mousse for Children



What is RECALDENT™ (CPP-ACP)?

Nature has developed highly efficient protein systems in saliva, blood and milk that stabilise and transport calcium and phosphate, which are the fundamental building blocks of our teeth and bones. Cow's milk is recognised as having the most effective system for delivery of calcium and phosphate. Extensive research over the past 25 years has refined the specific milk peptide which transports these essential minerals to help build, strengthen and protect our teeth – casein phosphopeptide amorphous calcium phosphate – otherwise known as RECALDENT™ (CPP-ACP).

Recaldent



Nature knows best

Made from milk



"The use of GC Tooth Mousse on deciduous teeth after they have erupted helps them mature and become more resistant to decay due to the large amounts of bio-available calcium and phosphate provided by the RECALDENT™ (CPP-ACP) contained in GC Tooth Mousse. GC Tooth Mousse also helps reverse the initial stages of decay, as sometimes our feeding habits and oral care of infants is not as good as we hoped, despite our best intentions. Also, there is no fluoride in the product, so swallowing small amounts is not a problem."

DR DAVID MANTON, UNIVERSITY OF MELBOURNE
CONVENOR OF PAEDIATRIC DENTISTRY... AND FATHER



This is the ideal time to reinforce good oral care habits started during infancy, which will have benefits for a lifetime. It is also a time when permanent teeth are erupting and additional protection, especially of erupting molars, may be required.

"Often erupting molars are hard to clean and are therefore susceptible to early decay, especially as the enamel has yet to mature. The use of GC Tooth Mousse daily can decrease this risk of decay whilst the tooth erupts, before it can be sealed effectively."

DR DAVID MANTON, UNIVERSITY OF MELBOURNE, CONVENOR OF
PAEDIATRIC DENTISTRY... AND FATHER

Mousse for the mum to be



"During pregnancy, hormonal changes result in a dramatic decrease in the ability of saliva to counteract acid attack. At the same time, nausea from morning sickness can lead to frequent reflux of acid contents of the stomach, particularly in the second and third months of pregnancy. Reduced protection from saliva, combined with frequent acid reflux, means that many pregnant women are at high risk of tooth decay and erosion during pregnancy. Regular use of GC Tooth Mousse during pregnancy helps to maintain elevated levels of calcium and phosphate in the saliva, protecting the teeth from dental erosion. GC Tooth Mousse will also help soothe sensitive teeth and protect areas of decay."

PROF LI WALSH, UNIVERSITY OF QUEENSLAND

Ortho Mousse



"GC Tooth Mousse can help prevent plaque accumulation around brackets, archwires, springs and other appliances, and can assist in countering the effects of acids produced by dental plaque."

Regular application of GC Tooth Mousse throughout the full course of orthodontic treatment can prevent areas of decalcification (white spots) developing."

PROF LI WALSH, UNIVERSITY OF QUEENSLAND

Application guide

Indication	How often	Duration
General prevention	Once a day, after flossing and brushing in the evening	Ongoing as part of a general prevention program
Additional comments: Regular use of GC Tooth Mousse will help maintain a healthy oral environment		
Active decay	Twice a day, after flossing and brushing	Until risk of future decay has been reduced
Additional comments: Your dental professional can assess and determine the potential source(s) of risk and can make suggestions on how to reduce risk		
White spot lesions	Twice a day, after flossing and brushing	8-12 weeks and thereafter as required
Additional comments: There are different types of white spots and your dental professional may need to pre-treat the white spot surface prior to starting application of GC Tooth Mousse		
Orthodontic treatment	Twice a day, after flossing and brushing	During entire orthodontic process
Additional comments: Complete treatment 12 weeks after the finish of orthodontic procedures, or thereafter as required to reverse any white spot lesions		
Developmental defects in enamel	Minimum twice a day, after flossing and brushing, and as required for sensitivity	Continuous treatment as required
Additional comments: Depending on severity, your dental professional may need to undertake additional procedures to protect these teeth		