

OMNILUX clear-U™ has been developed from a medical device that has been used by dermatologists worldwide for the successful treatment of acne.

What is clear-U?

OMNILUX clear-U™ uses light technology in the form of light emitting diodes to:

- Destroy the bacteria responsible for acne
- Stimulate healthy tissue growth
- Decrease inflammation and pore size
- Limit oil production in the future

When should I use clear-U?

Use OMNILUX clear-U™ at the first sign of an acne breakout. The treatment program consists of 8 treatment sessions spread over 4 weeks, 2 treatments per week. Use the blue light on the first treatment session; wait 2 days, then use the red light. Repeat this for the next 3 weeks.

Who can benefit from clear-U?

Clear-U benefits anyone suffering with mild to moderate inflammatory acne.

How does the treatment work?

One of the key bacteria responsible for acne is *Propionibacterium acnes* (P.acnes). This bacteria produces natural chemicals called porphyrins, which are sensitive to light at specific wavelengths. When your skin is flooded with blue and red light, the bacteria are killed and inflammation is decreased. Your acne fades and your skin heals. It's as simple as that!



ACTUAL SIZE: 5" h x 2.5" w x 1.5" d, 2.9oz.



Photos courtesy of Lee SJ, Park KH, Department of Dermatology, NMC, Korea

“The Clear-U handheld unit provides our patients with the tried and true results of LED technology in the palm of your hand. Finally an at home unit that really works to treat acne! LED technology is effective and easy to use. What more could you ask for?”

What results can be expected from a series of treatments?

Best results will be seen between 4-8 weeks after treatment. The light stimulates natural processes that occur after the treatment has stopped. Studies have shown complete clearance of lesions, however, on average you should expect 7 out of 10 acne lesions to clear.

How many treatments will I need?

The course of treatments is 2 x 20 minute treatments per week for 4 weeks of alternating blue and red light. However, because the light will not damage your skin, there is no limit on the number of treatments you can have.

How long do the effects last?

For most people, their acne is clear for six months. In some cases, follow-up treatments (typically 2 to 4 treatments) may be necessary.

Is clear-U all I need for my acne?

OMNILUX clear-U™ is not a substitute for good skin hygiene or a doctor's care. We suggest that the light therapy is used in conjunction with a preparation containing salicylic acid. Your skin specialist can discuss this combination treatment with you.