

# Motobility

*You've got a friend*



**CUSTOM MADE ADJUSTABLE BEDS**

[www.motobility.com.au](http://www.motobility.com.au)



# COME AND SAY HELLO!



## Osborne Park

📍 311 Selby Street North, Osborne Park, WA 6017

☎ 08 9242 7333



## Joondalup

📍 31/57 Joondalup Drive, Edgewater, WA 6027

☎ 08 9220 9040



## Mandurah

📍 2/17 Gordon Road, Mandurah, WA 6210

☎ 08 9479 3232



[www.motobility.com.au](http://www.motobility.com.au)



[info@motobility.com.au](mailto:info@motobility.com.au)



[www.facebook.com/motobility](http://www.facebook.com/motobility)



[www.youtube.com/@motobilityAU](http://www.youtube.com/@motobilityAU)



## CONTENTS

---

<b>01</b>	About Us	<b>02</b>	Our Custom Made Beds
<b>03</b>	Benefits of an Adjustable Bed	<b>04</b>	How to Choose an Adjustable Bed
<b>06</b>	Bed Functions	<b>10</b>	Peak Home Care Bed
<b>12</b>	Pressure Care and Bed Sores	<b>13</b>	Mattress Comparison Guide
<b>14</b>	Types of Mattresses	<b>22</b>	Bed Sizes and Measurements
<b>24</b>	Bed Surrounds	<b>26</b>	Pillows, Manchester and Accessories
<b>28</b>	Servicing and Government Grants	<b>29</b>	Other Equipment



Free in-home  
test drives



Australian Owned  
& Operated



Mobile Servicing  
& repair



Friendly &  
Knowledgeable Staff



Roadside  
Assistance



# IMPORTED AND WA CUSTOM MADE ADJUSTABLE BEDS



We have an extensive range of adjustable beds available and have the look and feel of a sophisticated ensemble but with all the benefits and ease of a hospital bed. We custom make all beds on-site here in Perth, at our Osborne Park premises, using high-grade carbon steel frames, made with higher density steel than many others on the market for superior strength and longevity.

The motors are manufactured in Germany to hospital standards. This means they are sealed against dust and water, and for your safety, protected by a thermal overload switch. Based on your personal needs, we can fit different motor options, including back-lift, leg-lift and massage. Motobility can cater for all regular bed sizes as well as custom sizes.



# **BENEFITS OF AN** ADJUSTABLE BED

## **FOR BACK PAIN**

When you sleep on adjustable beds, you will be able to do so in a position that more closely matches the contours of your body and this in turn will mean that you are able to keep your back flatter against the ground. This can help to prevent conditions such as sciatica. By supporting the base of the spine it is possible to keep the spine aligned and help to prevent trapped nerves and release those that are currently lodged.

At the same time, by having some incline you can reduce the pressure on your back that causes it to compress during the night. This is useful for preventing yourself from waking up and causing injury to yourself while the spine is more vulnerable. You'll also be able to support your head well using adjustable beds, and this will mean you don't need to pile up too many pillows which can cause neck ache and headache.

## **FOR APNEA AND SNORING**

If you snore then using an adjustable bed, it is possible to use the incline to help improve your health. This is because one of the major causes of snoring is that your windpipe closes over due to the weight of your neck and this then results in the air making a snoring noise as it passes through. By propping yourself upright slightly with adjustable beds, it is possible to alter the direction of the gravity so that it does not obstruct your windpipe. From there, you can then reduce the sound of snoring and improve the quality of your sleep. This is particularly important for those who snore due to sleep apnea.

## **FOR DIGESTION**

Being slightly upright is also highly beneficial for your digestion as it helps your body to process food during the night. It is not advisable to go to sleep on too full a stomach at any point, but if you have eaten within the last few hours then addressing your position can be helpful and by using adjustable beds you can position yourself at the recommended 6 inches and this will help to encourage healthy digestion.

## **FOR SWELLING**

If ever you sprain or pull a limb, break it, or otherwise hurt yourself, then you will be told to elevate it in order to prevent the pain. Often you will use a pillow or cushion to attempt this. However, you are likely to find in doing so that you move off during the night. On the other hand, by using adjustable beds, you can give yourself a better elevation and thereby help reduce swelling.

Those who have jobs working on their feet can significantly benefit from this if they often come home and find their lower legs and feet swollen. By elevating their legs slightly, it's possible to not only feel much better but also to help combat problems.

# HOW TO CHOOSE AN ADJUSTABLE BED

1



## FUNCTIONS

The first step in choosing an adjustable bed is to identify what functions you will need. This allows you to get the most out of your adjustable bed and ensure a good nights sleep.

- The head lift function helps people to breathe easily when their head/back is elevated at an angle.
- Dual lift adjustable beds are very sought after when people are looking for a good nights sleep and help alleviate the bodies aches and pains. With the dual lifts working together, the goal is to achieve the 'zero-gravity' position where your body feels most comfortable. It helps to redistribute the weight of your body evenly rather than focusing all on one pressure point. If you have swollen feet, it is ideal for elevating your feet slightly higher than your body which helps prevent swelling.
- The Hi-Lo function is beneficial for those that struggle to get in and out of bed. This is also beneficial for carers and nurses as it reduces the amount of bending over and strain on their back.

2



## MATTRESSES

Selecting a mattress can be a very personal selection depending on your specific needs. Some people may require a high-pressure care mattress due to existing medical conditions. Your medical practitioner or occupational therapist should always advise you before purchasing a mattress.

3



### SIZE

Choosing the right size is not too difficult as most people already know what they want or are simply replacing what they currently have at home. However, if you have a special need or a medical condition which might require a particular sized bed to suit. This should always be advised by your medical practitioner or occupational therapist to ensure choosing the correct sized bed for you.

4



### BED SURROUNDS

Headboards and surrounds enable you to customise the look and style of your bed. In the NDIS dictionary, the headboard and surrounds are deemed as non-essential items. Therefore, it is classified as an optional extra. Our custom made surrounds allow you to choose the colour, type of fabric and style.

5



### ACCESSORIES & MANCHESTER

Bed Accessories can make a world of difference on how you use your bed on a day to day basis. While you may not need any accessories, it's good to keep the selection in mind in case your needs change. Sometimes we don't realise how much a small aid can help assist us with our daily living. Motobility's specially trained team work together with the Health industry and Occupational Therapists to listen to your needs and advise you on which accessory product will work best for you.



# HEAD LIFT FUNCTION




  
REGULAR

  
HEAD LIFT

## FEATURES

- Will raise the head of the bed up to 60 degrees.
- Can be used in conjunction with the Hi-Lo function, making it a Headlift Hi-Lo.
- Lifting the head of the bed allows for easier entering/exiting of the bed.
- Perfect for reading, watching TV and sitting up to eat.
- Massage function available

Minimum Base Height (Floor to top)	250mm
Maximum Base Height (Floor to top)	No Limit
Minimum Bed Height (Including Mattress)	450mm
Maximum Bed Height (Including Mattress)	No Limit
Standard Bed Height	565mm
Maximum Lift	N/A
Leg Elevation	N/A
Head Lift Angle	60°



AVAILABLE IN ALL MATTRESS SIZES



# DUAL LIFT FUNCTION



REGULAR



HEAD LIFT



LEG LIFT




ZERO GRAVITY



## FEATURES

- Allows you to independently control the backrest of the bed and feet of the bed.
- Head-lift will raise the Head of the Bed up to 60 degrees.
- Leg-lift includes a knee break, lifting the legs from the knees for additional comfort and allow a raise of up to 270mm.
- Can be used with the Hi-Lo function, making this bed a Hi-Lo Dual Lift.
- Allows you to reach a maximum level of comfort in your bed.
- Massage function available

Minimum Base Height (Floor to top)	250mm
Maximum Base Height (Floor to top)	No Limit
Minimum Bed Height (Including Mattress)	450mm
Maximum Bed Height (Including Mattress)	No Limit
Standard Bed Height	565mm
Maximum Lift	N/A
Leg Elevation	0-270mm
Head Lift Angle	60°



AVAILABLE IN ALL  
MATTRESS SIZES

# HI-LO FUNCTION



  
REGULAR


  
HI-LO



## FEATURES

- Will lift the bed horizontally by up to 340mm.
- Can be used in conjunction with the Dual lift function and head-lift function.
- This function is required by some in-home care services due to Occupational Health and Safety of their staff.
- The horizontal lift function can be used as assistance in entering and exiting your bed. Suitability depends on the type of mattress and the abilities of the user.
- Massage function available

Minimum Base Height (Floor to top)	330mm
Maximum Base Height (Floor to top)	No Limit
Minimum Bed Height (Including Mattress)	530mm
Maximum Bed Height (Including Mattress)	No Limit
Standard Bed Height	565mm
Maximum Lift	340mm
Leg Elevation	N/A
Head Lift Angle	60°



**AVAILABLE IN ALL  
MATTRESS SIZES**

# SPLIT SYSTEM ADJUSTABLE BED




## FEATURES

- Perfect for couples.
- Available in Split-Queen and Split-King sizes.
- Allows each individual to choose a separate mattress and functions.
- Companion beds are available, which have no functions but can have functions retro-fitted at any given time in the future.
- Great option if one individual has higher care needs than the partner.
- Companion Beds are a great way to reduce costs.
- Massage function available



2 x Long Small Singles = Split Queen



2 x Long Singles = Split King

# PEAK HOME CARE BED



## AVAILABLE FUNCTIONS



REGULAR



HI-LO



HEAD LIFT



LEG LIFT



TRENDELENBURG & REV



## FEATURES

- Fully featured home care style electric bed.
- Durable powder coated frame.
- Australian Standards compliant.
- Linak Motor system with battery back-up.
- 10 Years Limited warranty (frame)
- 2 Years Full Warranty (Electrics, Motor)
- Full range of accessories available.
- Bed surrounds available in two fabric colours - Midnight and Mist.

**Available in King Single size**



## HANDSET

Simple 10 button  
backlit handset.





## ACCESSORIES

**High Side Bed Rail**



**Low Side Bed Rail**



**Clamp On Bed Stick**



Max bed height	73cm
Min bed height	26cm
Length	213cm
Width	117cm
Trendelenburg & Rev	-15° to +15°
Maximum Patient Weight	180kg
Product Weight	94kg
Recommended mattress length	203cm
Recommended mattress width	107cm
Backrest adjustment	0 to 70°
Knee-break adjustment	0 to 35°

## BED SURROUNDS COLOURS

**MIDNIGHT**  
(SJ5002)

**MIST**  
(SJ5001)

# PRESSURE CARE AND BED SORES

## WHAT ARE PRESSURE SORES?

These annoyingly uncomfortable sores (also known as pressure ulcers) are injuries resulting from prolonged pressure and result in injuries to the skin and underlying tissue. Our skin is much more delicate than we know, placing too much pressure in any given place for an extended period of time can cause damage. Most pressure sores tend to show up in bony areas of the body such as the hips, buttocks, tailbone, ankles and heels. As a result, the skin can become thin, flaky, puffy, clammy, discoloured or even broken. Unfortunately, these sores can develop quickly, and they can be tough to get rid of. For this reason, one of the best things to do is actually to prevent them from happening in the first place.

### 1. PRESSURE

Constant pressure on any part of your body can lessen the blood flow to tissues. Blood flow is essential to delivering oxygen and other nutrients to tissues. Without these essential nutrients, skin and nearby tissues are damaged and might eventually die.

### 2. FRICTION

Friction occurs when the skin rubs against clothing or bedding. It can make fragile skin more vulnerable to injury, especially if the skin is also moist.

### 3. SHEAR

Shear occurs when two surfaces move in the opposite direction. For example, when a bed is elevated at the head, you can slide down in bed. As the Tail-bone moves down, the skin over the bone might stay in place — essentially pulling in the opposite direction.

## PRESSURE MANAGEMENT

An essential part of bed sore prevention involves using the right mattress. Pressure mattresses are designed to spread the weight of your body over a larger surface area, offering a more comfortable night's sleep and a reduced risk for sores to develop. Pressure-relieving foam mattresses provide what's called a "reactive surface" and generally contain multiple layers of special foam that "react" or "respond" to a person's body, heat, shape and movement. Pressure can also be relieved using pockets of air created in the layers of the foam.



#### Low Risk


































If users are not prone to pressure sores, low risk ratings are suitable








#### Very High Risk

Ideal for users that have a high risk of getting pressure sores

# MATTRESS COMPARISON GUIDE

	Perfect Support	Luxury Support	Response Contour	Flexiblock	Quantam Comfort	Natures Contour
Main Material	 <p>Bonnell Spring</p>	 <p>Bonnell Spring with Pillowtop</p>	 <p>Pocket Spring/ Gel Infused Latex</p>	 <p>Memory Foam</p>	 <p>Gel Infused Memory Foam</p>	 <p>Gel Infused Latex</p>
Feel	 <p>Medium</p>	 <p>Medium Plush</p>	 <p>Med/Firm</p>	 <p>Firm or Medium</p>	 <p>Plush</p>	 <p>Plush</p>
Quilting	Wool Blend	Wool Blend	Wool Blend	Carrflex Waterproof Cover	Wool Blend	Wool Blend
Height	200mm	200mm	200mm	200mm	200mm	200mm
Edge Support	Soft	Medium	Strong	Strong	Medium	Medium
Recommended Sleeping Positions	 <p>Back Sleeper</p>	 <p>Back Sleeper</p>	 <p>Tummy Sleeper</p>	 <p>Tummy Sleeper</p>	 <p>Side Sleeper</p>	 <p>Active Sleeper</p>
	 <p>Active Sleeper</p>	 <p>Active Sleeper</p>	 <p>Back Sleeper</p>	 <p>Back Sleeper</p>	 <p>Min. Partner Disturb</p>	 <p>Side Sleeper</p>
			 <p>Min. Partner Disturb</p>	 <p>Min. Partner Disturb</p>		 <p>Min. Partner Disturb</p>
Pressure Rating	 <p>Low Risk</p>	 <p>Low Risk</p>	 <p>Low to Medium Risk</p>	 <p>Very High Risk</p>	 <p>High Risk</p>	 <p>Medium Risk</p>

# PERFECT SUPPORT MATTRESS


Main Material	 <p>Bonnell Spring</p>
Feel	 <p>Medium</p>
Quilting	Wool Blend
Height	200mm
Edge Support	Soft
Recommended Sleeping Positions	 <p>Back Sleeper</p>
	 <p>Active Sleeper</p>
Pressure Rating	 <p>Low Risk</p>

## PERFECT SUPPORT

Our Perfect Support is a Bonnell Orthoflex Innerspring mattress, with Ortho 3-zone support to assist with correct spinal alignment.

It offers a nice bounce and cloud like comfort feel, yet offers a lot of support all with a medium plush feel.






## MATTRESS INTERNALS



- 1 Wool Blend Comfort Firm Quilting
- 2 Foam Layer
- 3 Everlast Insulation Pad
- 4 Bonnell 3-Zone Orthoflex Springs
- 5 Everlast Insulation Pad



# LUXURY SUPPORT MATTRESS

Main Material	 Bonnell Spring with Pillowtop
Feel	 Medium Plush
Quilting	Wool Blend
Height	200mm
Edge Support	Medium
Recommended Sleeping Positions	 Back Sleeper
	 Active Sleeper
Pressure Rating	 Low Risk

## LUXURY SUPPORT







Our Luxury Support builds on the Perfect Support. Utilizing a Bonnell Orthoflex Luxury Support Innerspring system, we've added a comfort pillow-top as well as a Edge Care Comfort Plus foam box system, creating a luxurious comfort feel and maximizing the sleep surface and decreasing the risk of rolling out of beds. Slightly less bouncy, yet very supportive.

## MATTRESS INTERNALS



- 1 Wool Blend Comfort Firm Quilting
- 2 Body Contouring Comfort Plush
- 3 Body Contouring Enhanced Comfort
- 4 Edge Care Comfort Plus
- 5 Slimline Bonnell Orthoflex Springs
- 6 Foam Layer

# RESPONSE CONTOUR MATTRESS

Main Material	 Pocket Spring/ Gel Infused Latex
Feel	 Med/Firm
Quilting	Wool Blend
Height	200mm
Edge Support	Strong
Recommended Sleeping Positions	 Tummy Sleeper
	 Back Sleeper
	 Min. Partner Disturb
Pressure Rating	 Low to Medium Risk

## RESPONSE CONTOUR






This mattress is a hybrid of two of the most popular mattresses a spring system and latex. The pocket spring system has been engineered to intuitively respond to your bodyweight and shape, for individual support and minimal partner disturbance. The Natural Gel-infused Latex provides the comfort and health benefits of a Latex mattress.

## MATTRESS INTERNALS



- 1 Wool Blend Comfort Firm Quilting
- 2 Latex Gel Infused
- 3 Comfort Support Plus
- 4 Edge Care Comfort Plus
- 5 Contour Response Pocket Spring
- 6 Comfort Support Plus


# QUANTUM COMFORT MATTRESS

Main Material	 Memory Foam
Feel	 Plush
Quilting	Wool Blend
Height	200mm
Edge Support	Medium
Recommended Sleeping Positions	 Side Sleeper
	 Min. Partner Disturb
Pressure Rating	 High Risk

## QUANTUM COMFORT







Similar to our Flexiblock mattress, yet combined with our gel-infused memory foam and available with a regular quilting. The Quantum Comfort gives the utmost in satisfaction, both for your health and comfort. You don't lie on this mattress... It holds you. Also infused with Fusion Gel, the Quantum Comfort will assist with helping to regulate body temperature.

## MATTRESS INTERNALS



- 1 Wool Blend Comfort Firm Quilting
- 2 Fusion Gell Layer
- 3 Edge Care Comfort Plus
- 4 Body Contouring Enhanced Comfort
- 5 Quantum Comfort Layer

# NATURES CONTOUR MATTRESS

Main Material	<div>  <p>Gel Infused Latex</p> </div>
Feel	<div>  <p>Plush</p> </div>
Quilting	Wool Blend
Height	200mm
Edge Support	Medium
Recommended Sleeping Positions	<div>  <p>Active Sleeper</p> </div>
	<div>  <p>Side Sleeper</p> </div>
	<div>  <p>Min. Partner Disturb</p> </div>
Pressure Rating	<div>  <p>Medium Risk</p> </div>

## NATURES CONTOUR

The Natures Contour Mattress is the Ultimate adjustable bed mattress. It conforms and supports the body ergonomically without resisting body weight instead welcomes it and makes you feel like you are lying on the mattress, rather than in it. Our natural latex has been infused with gel particles to increase it’s temperature regulation properties.







## MATTRESS INTERNALS



- 1 Wool Blend Comfort Firm Quilting
- 2 Body Contouring Enhanced Comfort
- 3 Latex Gel Infused
- 4 Foam Layer




# FLEXIBLOCK MATTRESS

Main Material	 Memory Foam
Feel	 Firm or Medium
Quilting	Carrflex Waterproof Cover
Height	200mm
Edge Support	Strong
Recommended Sleeping Positions	 Tummy Sleeper
	 Back Sleeper
	 Min. Partner Disturb
Pressure Rating	 Very High Risk

## FLEXIBLOCK






The Flexiblock mattress is an entry level pressure relieving mattress. It comes with a waterproof cover, which offers a breathable surface that also assists in minimizing shearing. It is built out of 4 different types of foam and has a 50mm memory foam topper. The waterproof cover is made with a toxic-free, safe and environmentally friendly material called polyurethane. The polyurethane material is also 100% latex-free, which is ideal for those with latex allergies. The waterproof membrane protects against body fluids that can cause a damp mattress, mould and fungi, which can all contribute to a toxic sleeping environment.

## MATTRESS INTERNALS



- 1 Waterproof Membrane
- 2 Pressure Relieving Foam
- 3 Plush Comfort Foam
- 4 Edge Care Comfort Plus
- 5 Weight Absorbing Foam


# METZELLER PRESSURE CARE MATTRESS

Material	Made from eco-friendly materials like sunflower oil, instead of petroleum products
Feel	 Soft or Medium
Quilting	Waterproof Co
Height	180mm
Edge Support	Strong
Recommended Sleeping Positions	<div>   Side Sleeper         </div> <div>   Back Sleep         </div> <div>   Min. Partner Disturb         </div>
Pressure Rating	 Very High Risk

## METZELER PRESSURE CARE

Incorporating patented Rubex™ Duo Cell technology, a Metzeler Mattress will provide you with the perfect night’s sleep. A Metzeler mattress offers superior comfort and support which provides optimal pressure distribution, ventilation and spinal alignment.


- Optimal temperature management
- Multi-zone support system
- Durable reinforced mattress edge
- Waterproof breathable cover
- Environmentally friendly materials



- The only static mattress to offer vertical and horizontal pressure support
- Made of natural non-toxic environmentally friendly raw materials
- Optimises ventilation, breathability and enhances thermoregulation
- Offers optimal comfort and support for individuals with a disability, including those with Central Nervous System (CNS) disorders, post-surgery or cancer treatment, and conditions that may impact thermoregulation such as CP, MND and MS,
- Can be used on both sides



## MATTRESS INTERNALS



- ① **Waterproof, breathable cover that's moisture-managing, easily removable, and washable at 95°C**
- ② **Mattress structure and materials provides superior thermoregulation**
- ③ **Multi-zone system offers optimal postural support and pressure distribution**
- ④ **Reinforced foam edges**
- ⑤ **Weight Absorbing Foam**

# BED SIZES AND MEASUREMENTS



**Standard Single**  
920mm x 1880mm  
3' x 6'1" (inch)



**Standard Double**  
1370mm x 1880mm  
4'5" x 6'1" (inch)



**Long Small Single**  
762mm x 2030mm  
2'6" x 6'8" (inch)  
x2 = Split Queen



**Long Single**  
915mm x 2030mm  
3' x 6'8" (inch)  
x2 = Split King



All mattresses have a 200kg max weight capacity





### King Single

1067mm x 2030mm  
3'6" x 6'8" (inch)



### Long Double

1371mm x 2030mm  
4'6" x 6'8" (inch)



### Queen

1524mm x 2030mm  
5' x 6'8" (inch)



### King

1824mm x 2030mm  
6' x 6'8" (inch)



All mattresses have a 200kg max weight capacity

# BED SURROUNDS



Whilst all our Beds are freestanding and have four lockable heavy duty casters for easy movement, you may want to choose a bed surround to match your current interior. We offer a range of upholstered bed surrounds, including WA custom-made designs, which will give you a lot of creative freedom to achieve your desired look. We also offer custom-made head and foot-boards.





## MOUNTED HEADBOARD AND FOOTBOARD

The Headboard and foot-board attach directly to the bed. This allows you to easily move the bed to a different location or room. If you have the Hi-Lo function fitted the headboard and foot-board will move up and down with the bed.



## BED SURROUNDS AND HEADBOARDS

We offer a range of different upholstered bed surrounds, including budget friendly stocked items and WA-made custom options. All our bed surrounds are free-standing and aren't attached to the bed. They can easily be dismantled thanks to the clever click in place design. Choose from box style surrounds and open sides for more under bed access.



## NO SURROUND (THE NAKED BED)

This is the standard finish of the bed. It shows the black metal frame, the grey upholstered bed base, which adjusts the mattress and the mattress that attaches via a horizontal velcro strip. This style can easily be dressed up, by using an oversized dooner or quilt.



## IMPORTED SURROUNDS

Our imported surrounds come in beige or dark grey.



# PILLOWS AND MANCHESTER

Many people are not aware of what a significant role selecting the right pillow plays. The correct pillow should provide support to your neck and head and shouldn't "deflate" or get lumpy. Having the wrong pillow can cause tension headaches, neck pain breathing issues, such as snoring, acid reflux and even allergy related issues due to dust mites.

We provide a large range of solutions to provide a moisture barrier and protect the mattress from excessive moisture as well as incontinence. Our Staydry range offer great solutions for everyday use that can easily be laundered and offer a natural, breathable feel.



**Terry Towelling  
Mattress Protector**



**Zip-on Waterproof Mem-  
brane (Dartex Cover)**



**Baronne Quilted Mattress  
Protector Fitted Style**



**Sensi Select Bamboo Charcoal  
Infused Memory Foam Pillow**



**King Koil Pocket Spring  
Pillow – Low Profile**




**Sensicloud Rebound Smart  
Foam Pillow**



**King Koil Pocket Spring  
Pillow – High Profile**



**Sensi Polar Active Touch Gel  
Infused Memory Foam Pillow**



**Sensi Gel Infused Memory  
Foam Pillow**



# BED ACCESSORIES

Bed Accessories can make a world of difference on how you use your bed on a day to day basis. While you may not need any accessories, it's good to keep the selection in mind, in case your needs change. Sometimes we don't realise how much a small aid can help assist us with our daily living. Motobility's specially trained team work together with the Health industry and Occupational Therapists to listen to your needs and advise you on which accessory product will work best for you.



**Bed Rail Hoop Design**



**Under Mattress Bed Pole**



**Under Mattress Bed Rail**



**Collapsible Side Rail**



**Over Bed Table**



**Bed Cradle**



**Mattress Holding Bar**



**Battery Backup**



**Angled Heavy Duty Bed Rail**

# SERVICING



In order to uphold the warranty\* an annual service of the bed is required. This service is to be carried out by one of our friendly and well trained technicians.

\*Refer to the product's Warranty & Maintenance Guide for warranty information.

## GOVERNMENT GRANTS

### WHAT ARE GOVERNMENT GRANTS?

Over the past few years, there has been a big push by the government to promote “ageing in place”, which means trying to keep people in their own homes or in retirement communities, rather than moving to a nursing home. In order to achieve this, the government is offering support, including services and equipment grants. Depending on your circumstances, you may be eligible.

### WHAT HELP CAN I RECEIVE?

With the introduction of the [NDIS \(National Disability Insurance Scheme\)](#), the sector has undergone the most significant reform since Medicare. While there have been mixed reviews, the NDIS promises a more longterm approach of assistance and education. It is important to note that you have to be 65 years or younger to be considered for the NDIS. If you are over 65 years of age, you may be eligible for one of the below.

[Home Care Packages \(HCP\)](#) through My Aged Care are a way that older Australian's can access affordable services to get assistance around the home. These are designed to cater for more complex needs and can include allowances for Assistive Technology (AT) at higher levels.

[Equipment For Living Grant \(EFL\)](#) through the Independent Living Centre ILC) funded by the Department of Communities (formerly DSC) aims to assist people in obtaining vital equipment not currently available through other funding channels. These need to be essential items that will support the user in everyday tasks around their homes and the community.

# OTHER EQUIPMENT AVAILABLE FROM MOTOBILITY



## SCOOTERS & ELECTRIC WHEELCHAIRS

Motobility offer an extensive range of travel scooters, portable scooters and electric wheelchairs from the mobility industry's most reputable brands including Solax, Luggie, Supa Scoota, E-Traveller and more. Choose from lightweight automatic folding travel scooters, electric wheelchairs, manual-folding scooters, large scooters and portable pull apart scooters.



## PLATFORM LIFTS

Stairs should not stop anyone from entering a public building, a place of employment or stop them from continuing to live in their own home. Motobility is the exclusive Australian distributor of wheelchair platform lift market leader Lehner Liftechnik. An Austrian company, Lehner manufactures two types of wheelchair platform lifts which both offer a safe and reliable solution to climbing stairs. The simple design allows easy operation from the handset or the remote control.



## STAIRLIFTS

A staircase should not stop anyone from continuing to live in their own home. Motobility is the local distributor of world-renowned Stannah Stairlifts in Australia; the safe, reliable and painless alternative to climbing the stairs. Stannah applies a simple design that provides easy operation from the seat mounted joystick or remote control. All stair lifts are fully installed and maintained by our trained technicians.



## LIFT AND RECLINE CHAIRS

Motobility's individually handcrafted Lift & Recline Chairs are made according to your seat size, leg length and the shape of your back with an ergonomic design featuring a full leg chaise to give you the utmost comfort while in a seated or reclined position, and a reliable motor assisting you as you rise from it.





# Motobility

*You've got a friend*

[www.motobility.com.au](http://www.motobility.com.au)



SCAN TO SUBSCRIBE TO OUR  
**YOUTUBE**



## Joondalup

📍 31/57 Joondalup Drive, Edgewater, WA 6027

☎ 08 9220 9040

## Osborne Park

📍 311 Selby Street North, Osborne Park, WA 6017

☎ 08 9242 7333

## Mandurah

📍 2/17 Gordon Road, Mandurah, WA 6210

☎ 08 9479 3232