


OPULENCE BY INTRACEUTICALS



OPULENCE

An Intraceuticals Opulence Infusion® is a unique experience providing instant results you can see. Skin looks brighter and more radiant with uneven skin pigmentation improved.

 **intraceuticals**
hyperbaric oxygen technology



BEFORE

intraceuticals²

An Intraceuticals Opulence Infusion® targets the problematic areas of your face with a non invasive procedure.

AFTER



Graphic dramatisation
- not intended to reflect
an actual case study.


An Opulence Infusion is soothing, relaxing

and suitable for all skin types.

In just one soothing treatment your skin
will look younger, brighter, more luminous.

Uneven skin tone is balanced,
skin looks and feels hydrated, rejuvenated,
radiant and fresh.

Skin texture after an Intraceuticals Opulence Infusion® treatment.



intraceuticals
hyperbaric oxygen technology

HOW THE INFUSION WORKS

An Intraceuticals Opulence Infusion® harnesses three elements to achieve instant and dramatic Results:

A. Oxygen under light pressure aids delivery of age defying actives to the deeper layers of the skin at the same time as cooling and calming the skin.

B. Unique Opulence Serum contains lightweight Hyaluronic Acid, a powerful combination of Vitamins and antioxidants, and a blend of botanical brighteners. Hyaluronic Acid helps to instantly hydrate the skin resulting in increased firmness, reduced fine lines, and improved facial contours. Botanical brighteners help balance skin tone and Vitamin C protects and defends against environmental damage and premature skin ageing.




C. Treatment Support Products provide essential layers to seal in and supplement skin nutrients for maximum synergistic effect and longer lasting results.

HYALURONIC LAYERING

Lightweight Hyaluronic Acid in the Opulence Serum is absorbed quickly for maximum dispersion of nutrients. Opulence Hydration Gel with multi-weight Hyaluronic Acid offers a reservoir of nutrients for the lower levels of the skin, and Opulence Moisture Brightening Cream with multi-weight Hyaluronic Acid moisturises the skin's surface. The ultimate in intense hydration is now achieved and we call this Hyaluronic Layering.

Intraceuticals 3 Step Hyaluronic Layering process.



STEP 1: THE OPULENCE INFUSION

This is the first step in your skin's transformation. Cooling, calming, pressurised oxygen increases absorption of the Opulence serum and naturally improves your skin's moisture levels while vitamins and antioxidants protect and defend. A combination of lightweight Hyaluronic Acid, Vitamins A, C, E and Green Tea result instantly in smoother skin with reduced fine lines and wrinkles, better volume and increased firmness. A blend

of botanical brighteners and super concentrated Vitamin C help control melanin formation and balance uneven skin tone. Your skin will emerge deeply hydrated, more luminous, and glowing with health and vitality. Results are not only instant but continue to improve in the following days.



STEP 2: OPULENCE HYDRATION GEL

Intensify and maintain your results with Treatment Support Products. Opulence Hydration Gel with multi-weight Hyaluronic Acid provides the next hydration layer to seal in moisture and create an invisible barrier against moisture loss. Contains a blend of botanical brighteners including Alpha Arbutin and super concentrated Vitamin C to target uneven skin tone. Apply twice daily. Use Rejuvenate Contour Eye Gel at this stage too – specifically formulated for the delicate eye area with double the concentration of multi-weight Hyaluronic Acid.

STEP 3: OPULENCE MOISTURE BRIGHTENING CREAM

The final step in the hyaluronic layering process. This moisture rich product contains multi-weight Hyaluronic Acid providing intense hydration and a blend of botanical brighteners including Alpha Arbutin and super concentrated Vitamin C – to lock in moisture and restore skin's natural radiance. This Hyaluronic Layering process offers unmatched hydration and long lasting results. Continued use of Treatment Support Products is an essential stage in prolonging the benefits. Apply twice daily and prepare to be impressed.

OPULENCE TREATMENT SUPPORT PRODUCTS



Opulence Brightening Cleanser – 50ml (1.69 fl.oz.)

A gently foaming cleanser to thoroughly cleanse, lift away impurities and leave skin soft, clean and bright.
Use morning and night.

Opulence Hydration Gel – 40ml (1.35 fl.oz.)

A lightweight, non greasy, moisture rich gel for all skin types. Actively counteracts skin dehydration, the second step in Hyaluronic Layering and creates an invisible barrier against moisture loss.
Apply one to two pumps to cleansed skin morning and night.

Opulence Moisture Brightening Cream – 40ml (1.35 fl.oz.)

A light, nutrient rich cream and the third step in Hyaluronic Layering. Perfect for all skin types this soothing and healing nourisher locks in moisture and helps balance uneven skin tone.
Use after Opulence Hydration Gel morning and night.

Opulence Brightening Wand

A concentrated brightening complex to specifically target uneven pigmentation in a portable, take anywhere, anytime pen. Continued and regular application will help maintain luminous and even skin tone.
Apply twice daily or during the day as required.